


Product category: Sterile Dressing sets

REF: 610459

Description: Dressing pack no. 10



Technical drawing:



Material:

Code	Designation	Qty	Composition
5840	Forceps Burger 12,5 cm green	1	Polyethylene BG70 30% fiberglass
56005PC5	Swab gauze 8 ply 17 thr. 7,5 x 7,5 cm counted per 5	5	Absorbent gauze 100% cotton
56211	Ball gauze 17 thr. 20 x 20 cm Ø 30 mm	5	Absorbent gauze 100% cotton
1102	Drape PE/Tissue 40 x 45 cm blue	1	PE/Tissue 40 g/m2 (polyethylene 17 g/m2 + tissue paper 23 g/m2)
5844	Forceps Burger 12,5 cm blue	1	Polyethylene BG70 30% fiberglass

Package

Weight: 37 grams

Dimensions: 165 x 125 x 25 mm

Sterilization: Ethylene Oxide

Brand: Sterisets®

Packaging: Sterile single pack

Quantity in pcs	
Inner box	Transport carton
NA	44



Sterisets Medical Products
Zona Industrial 1 - Lote 11 a 14
4560-164 Guilhufe, Penafiel
Portugal www.sterisets.com

Conforms to ISO 13485
Conforms to Directive
93/42/EEC

Effective: 15-10-2019

Issued: DQ

Approved: DQ

REF 610459

MDDTF_N-5_11

ET 476/1

IMP 145/3